

## Soft Agar Assay Protocol

1. Preparation of Base Agar:
  - a. Dissolve 1% agarose (Difco Agar Noble) in sterile H<sub>2</sub>O, cool to 42°C in water bath. Use autoclaved 125 mL screwtop bottle
  - b. Warm 2X DMEM w/20% FBS and antibiotics to 42° in water bath.
    - i. It's best to equilibrate a 50 mL tube of DMEM and a 125 mL bottle of soft agar together in a water-filled 500 mL glass beaker within the water bath, so that they fit snugly and can be transported into the hood and kept at 37°C longer.
    - ii. 2X DMEM recipe: 38 mL 2X straight DMEM, 10 mL FBS (20%), 1 mL L-Glut (2X), 1 mL antibiotic (2X). I just used concentrations of FBS, L-Glut, and antibiotic that were already made up.
  - c. Equilibrate for about 45 minutes.
  - d. Bring beaker with DMEM and agarose into the hood.
  - e. Mix equal volumes of the agarose and the 2X DMEM to give 0.5% agarose/1X media. Also add appropriate amount of G418 if necessary.
  - f. Add 1 mL of mixture from step C to 6-well dish (1 mL/dish), tap sides of plate to spread evenly, allow 5 minutes to solidify.
    - i. Plates can be stored in 4° C for up to one week.
2. Preparation of Top Agarose:
  - a. Bring out Base Agar plates to equilibrate to room temperature if not already.
  - b. Dissolve 0.7% agarose in sterile H<sub>2</sub>O), cool to 37°C in water bath with 2X DMEM in the same manner as the base layer.
  - c. Count cells
    - i. Concentrations of cells I used were 2500, 7500, 20,000.
  - d. Add ½ media, ½ agarose solution. Make sure to add cells and G418 before adding the agarose, since the agarose quickly solidifies.
  - e. Add 1 mL to each well of the base agar layer, allow 15 minutes to solidify.
  - f. Add 1 mL 1X full media + G418 to each well.
  - g. Incubate plates at 37°C in humidified incubator 10-30 days (I did this over 14 days).
  - h. Feed cells twice per week.
3. Staining:
  - a. Stain cells with 0.5 mL Crystal Violet + methanol for >1 hour.
  - b. Wash Crystal violet off by adding 1-2 mL water to soft agar 4-5 times.
  - c. Image on a dissecting microscope.
  - d. Count colonies > 500 µm or other appropriate diameter. Easiest to use ImageJ → Binary → Analyze Particles if there are a lot of colonies.