

# Denaturing Purification His-Tagged Proteins From Mammalian Cells

1. Transfect cells with His-tagged construct
2. 36 hrs later, wash cells once in ice-cold 1X PBS
3. Collect cells in ice-cold 1X PBS
4. Add sample buffer to 10% of the cells to be used as input control
5. Remaining cells spin down at 4C
6. Resuspend in 6 mL Denaturing Lysis Buffer (containing BME)
7. Add 75 uL Ni-NTA Agarose beads
8. Incubate rotating at room temperature for 4 hrs
9. Wash beads for 5 min at room temperature successively in 750 uL each:
  - a. Denaturing Lysis Buffer without imidazole
  - b. Buffer A pH 8.0 without Triton X-100
  - c. Buffer A pH 6.3 + 0.2% Triton X-100
  - d. Buffer A pH 6.3 + 0.1% Triton X-100
10. Elute proteins with 75 uL Elution Buffer at room temperature for 20 min
11. Add 4X sample buffer to eluted proteins
12. Boil samples for 5 min
13. Run Input and Eluted protein samples on SDS-PAGE
14. Western Blot with appropriate antibodies

## **Denaturing Lysis Buffer**

6 M Guanidinium-HCl  
100 mM Na<sub>2</sub>HPO<sub>4</sub>-NaH<sub>2</sub>PO<sub>4</sub>  
10 mM Tris-HCl pH 8.0  
5 mM Imidazole (leave out in wash buffer)  
10 mM Beta-Mercaptoethanol (add fresh at time of use)

## **Buffer A**

8 M Urea  
100 mM Na<sub>2</sub>HPO<sub>4</sub>-NaH<sub>2</sub>PO<sub>4</sub>  
10 mM Tris-HCl pH 8.0 or 6.3  
10 mM Beta-Mercaptoethanol (add fresh at time of use)  
0.2% or 0.1% Triton X-100

## **Elution Buffer**

200 mM Imidazole  
150 mM Tris-HCl pH 6.7  
30% Glycerol  
5% SDS  
720 mM Beta-Mercaptoethanol (add fresh at time of use)