

His-Protein and IVT Binding Assay

- 1) Use 15 ug his-proteinX (test) and 15 ug of his-proteinY (neg. control) or no protein (neg. control)
- 2) Wash 14 uL of 50% Ni⁺⁺-NTA Qiagen beads in 400uL TBST, repeat (in 0.5mL tubes)
- 3) Add 50 uL of TBST (supplemented with 10mM B-Me) + His-tagged proteins of interest
- 4) Mix for 1 h in vibrating mixer (level 3-4) at RT (do binding in cold if protein is unstable at RT)
- 5) Remove supernatant
- 6) Add 50 uL of TBST + 5% milk + 10mM B-Me
- 7) Mix for 1 h in vibrating mixer (level 3-4) at RT
- 8) Remove supernatant (beads are hard to see, so be careful)
- 9) Add 50 uL of TBST + 5% milk + 10mM B-Me
- 10) Add 5 uL in vitro translated protein of interest. Keep 1 uL as input.
- 11) Mix for 1 h in vibrating mixer (level 3-4) at RT
- 12) Wash beads (red) 4X with 400uL TBST + 10mM B-Me
- 13) Remove all extra sup with pipette and add 20ul of 2X SDS-SB and boil 5-10 min
- 14) Also add 20ul 1X SDS-SB to input
- 15) Load 10 uL of input (10%), and 18 uL of bound
- 16) Run through stacker at 80 V and then 100 V until dye front reaches bottom of gel
- 17) Carefully remove gel and cut off stacking gel
- 18) Stain gel in coomassie blue stain at RT on rocker for 20 min or until gel is blue
- 19) Remove staining solution (reuse) and wash gel in dH₂O once
- 20) Destain gel in destaining solution with Kimwipes on rocker at RT for ~20 min or until bands are visible
- 21) Rinse gel in dH₂O
- 22) Dry gel on dH₂O soaked filter paper for 1 hr in gel dryer
- 23) Tape dried gel/filter paper down in phosphoimager cassette
- 24) Mark molecular weight ladder with excess ³⁵S-IVT reaction mixture
- 25) Place phosphoimager screen over gel and close cassette
- 26) Expose overnight
- 27) Image next day on Fuji PhosphoImager