

Ligation-Mediated PCR (LM-PCR)

Repair of Overhangs

1. Isolate genomic DNA using Qiagen DNAeasy spin column.
2. Use 10µg DNA per reaction. Add 3.3µL T4 DNA polymerase, 5µL 10X T4 polymerase buffer, 0.5µL 10mM dNTP (100µM final), and 0.5µL BSA. Make reaction up to 50µL with dH₂O.
3. Incubate for 15 min at 12°C.
4. Stop reaction with 1µL 0.5mM EDTA at 75°C for 20 min.
5. Purify DNA over Qiagen gel purification column. Elute in 30µL.

Linker preparation

1. Dilute oligos (BW-1 and BW-2) to 100 µM in dH₂O.
2. Use 2.5µL of each oligo in total of 50µL TENS buffer (250mM NaCl, 10mM Tris pH 8.0, 0.2mM EDTA).
3. Heat to 94°C for 2min, 70°C for 5min, 55°C for 5min, 50°C 5min, 45°C 5min, 40°C 5min, 35°C 5min, 30°C 5min, 23°C 5min, and finally on ice for 30 min. Always keep annealed oligos on wet ice, never greater temperature.
4. Dilute annealed linkers with 200µL TENS buffer for 1µM stock. Use 1µL of 1µM annealed linker per ligation reaction.

Ligation

1. Add: 10µL repaired genomic DNA, 1µL of 1µM annealed linker, 1µL T4 DNA ligase (4U/µL), 2µL 10X T4 ligase buffer, and 6µL dH₂O.
2. Incubate 16hr at 16°C in PCR machine.

Quantitative PCR (QPCR)

1. Perform QPCR. Use 2-5µL ligated, repaired genomic DNA with linker-2 and GFP#3 primers. Use GFP #7 primers as a normalization control for the amount of DNA in reactions.