

Making Competent Cells

Day 1:

1. Grow 5 mL LB culture of cells overnight at 37 °C
For XL1-Blue add 12.5 µg/mL Tetracycline

Day 2:

1. Dilute overnight culture into LB
For XL1-Blue – 1 mL of culture into 500 mL LB + 12.5 µg/mL Tet + 15 mM MgCl₂
For DH5α – 2 mL of culture into 120 mL LB
2. Grow cells at 37 °C until OD₆₀₀ = 0.4
3. Pellet cells at 3000 rpm for 15 min at 4 °C
4. Resuspend pellet in 150 mL cold MES A
5. Incubate on ice 20 min
6. Spin down at 3000 rpm for 15 min at 4 °C
7. Resuspend pellet in 30 mL of cold MES A + 15 % Glycerol
8. Aliquot (200-300 µL) and snap freeze in liquid nitrogen or dry ice ethanol bath
9. Transform cells – 100 µL for ligation or 50 µL for plasmid:
 - a. Add DNA – 10 µL of 20 µL ligation or 0.5 – 1.0 µL plasmid
 - b. Incubate on ice for 2 min
 - c. Heat shock 2 min at 42 °C
 - d. Incubate on ice for 2 min
 - e. Add prewarmed LB - 100 µL for ligation or 950 µL for plasmid
 - f. Plate - All 200 µL for ligation or 50 µL for plasmid

MES A (500 mL)

100 mL - 50 mM MES Salts pH 6.3 with KOH
2.77 g – CaCl₂·2H₂O
5 mL – 1M MnCl₂·4H₂O
dH₂O to 500 mL

Make fresh, filter sterilize, wrap in foil to protect from light, and store at 4 °C.