

Large Scale Immunoprecipitation

1. Harvest cells from 40 x 150 mm plates. Snap freeze and store in -80 °C
2. Lyse cells in 30 mL of ice-cold lysis buffer (50 mM Tris, pH 7.5, 0.2 M NaCl, 0.5% NP-40, 1 mM DTT, 1X PIC, 10 mM NaF, 5 mM β -glycerophosphate, 10 μ M cytochalasin B, 10% glycerol). Incubate on ice 15 min
3. Sonicate for 4 x 30 sec on ice and spin at 28,000 rpm with SW-28 for 45 min.
4. Filter lysate with 0.45 μ m filter
5. Split lysate into two tubes
6. Incubate half supernatant with 100 μ L anti-GST beads for 2 hr at 4 °C
7. Combine anti-GST depleted supernatant with other half of supernatant. Split into half again and add 100 μ L anti-MMS21 or 100 μ L anti-NSE1. Incubate 2 hr at 4 °C
8. Wash beads with 30 mL of lysis buffer once, then three times with 30 mL of wash buffer (50 mM Tris, pH 7.5, 0.4 M NaCl, 0.5% NP-40, 5 mM NaF, 2.5 mM β -glycerophosphate, 1 mM DTT), and finally once in 1 mL of lysis buffer
9. Transfer beads from last step in 1 mL of lysis buffer to eppendorf and elute with 150 μ L of glycine buffer (0.1 M glycine, pH 2.5, 1 mM DTT, 1X PIC) three times
10. Spin the elute to remove the beads and neutralize with 50 μ L of 1M Tris, pH 8.0.
11. Concentrate all 500 μ L to ~20 μ L. Remove supernatant.
12. Add hot SDS loading buffer to membrane and spin out
13. Run on one lane of 4-20% gradient SDS-PAGE gel
14. Wash gel three times with MiliQ water for 5 min each
15. Stain with Colloidal Blue overnight at 4 °C in gel plastic tray