

## Ubiquitin E2 Charging Assay

1. Setup thioester reactions in 1.5 mL tubes as follows:
  - a. 4  $\mu\text{L}$  of E1 (0.5  $\mu\text{g}/\mu\text{L}$ ) or 4  $\mu\text{L}$   $\text{dH}_2\text{O}$  as a control
  - b. 2  $\mu\text{L}$  His-UbcH2 (1  $\mu\text{g}/\mu\text{L}$ ) or 2  $\mu\text{L}$  UbcH2  $\Delta\text{Tail}$
  - c. 2  $\mu\text{L}$  10X thioester buffer
  - d. 6.5  $\mu\text{L}$  Flag-Ubiquitin (1  $\mu\text{g}/\mu\text{L}$ )
  - e. 5.5  $\mu\text{L}$   $\text{dH}_2\text{O}$
2. Incubate at room temperature for 10 min
3. Add 20  $\mu\text{L}$  2X non-reducing SDS sample buffer (homemade)
4. For reduced sample:
  - a. Remove 15  $\mu\text{L}$  into a new tube
  - b. Add 0.8  $\mu\text{L}$  of 1 M DTT (reduced sample only)
  - c. Boil for 5 min (reduced sample only)
5. For Non-reduced sample:
  - a. Incubate 15 min at 30  $^\circ\text{C}$
6. Run 15  $\mu\text{L}$  on SDS-PAGE gels. *Run at high voltage to prevent spontaneous reduction. Do not run reducing and non-reducing samples near each other, at least several lanes apart or on different gels.*
7. Coomassie stain gel

### 10X Thioester buffer

100 mM HEPES pH 7.5

1 M NaCl

400  $\mu\text{M}$  ATP

20 mM  $\text{MgCl}_2$